

# goldtrace

## The gold-plated fetal scalp electrode

### A message from the heart

Neoventa's Goldtrace® Fetal Spiral Electrode (FSE) is used when fetal heart rate monitoring during labour requires continuous and reliable internal registration. Goldtrace FSE is made of materials whose properties ensure the best-possible-quality signal for CTG.

Goldtrace FSE comes in a flexible protective guide tube that keeps the spiral tip inside the tube during application, protecting the user, mother and baby. The flexibility of the guide tube makes it possible to form according to desired anatomical curvature.

### Specifications

- Gold-plated spiral electrode
- Single use only
- Includes disposable applicator
- Not made with natural rubber latex
- Sterile product



### Article number

CNS000004/75 Goldtrace Fetal Spiral Electrode in boxes of 75 single packed pcs

### Materials

Spiral needle	316 Stainless Steel with 10-25µm Gold
Reference plate	316 Stainless Steel with 10-25µm Gold
Hub	PE
Connector	PE and brass, gold-plated
Wires	Copper with PVC jacket
Drive tube	PE
Guide tube	PE
Drive tube grip	PE

Follow the link for more information:  
<https://www.neoventa.com/products/goldtrace/>



PRM 400008/44 R5A © 2023 Neoventa Medical AB. All rights reserved.

Neoventa Medical AB  
Norra Ågatan 32, 431 35 Mölndal, Sweden  
+46 31 758 32 00, info@neoventa.com

[www.neoventa.com](http://www.neoventa.com)

neoventa   
A better start in life