







Part #	Identifier	Test Parameter	Intermediate cable	Compatible Device
4249		Tester for devices containing Masimo SET® or Masimo rainbow SET. Connects to all RD SET & Rainbow SET patient cables. Testers for SpO2 and PR only	4257; 4073; 4104; 4076; 4081	RAD 97, Radical 7, RAD 8, RAD 5 and OEM
4275		Tester for devices containing Masimo SET or Masimo rainbow SET, M20 instrument connector. Tests SpO2 , PR , SpCO , SpMet and SpHb	No cable	RAD 97, (RAD 67)
2368		Tester for devices containing Masimo SET® or Masimo rainbow SET. Tests SpHb , SpCO , SpMet , PVi , Pi , SpO2 and BPM .	No cable	Radical 7, RAD 8, RAD 5/5V
3502		Tester for devices containing Masimo SET or Masimo rainbow SET, 25-pin instrument connector. Tests SpO2 , PR , SpCO , SpMet and SpHb	Philips AO5 cables	Philips AO5 monitors
2367		9 pin tester. Tester for devices containing Masimo SET and connects to LNC patient cables. Tests SpO2 and PR only.	No cable	RAD 97, Radical 7, RAD 8, RAD 5 and OEM
1593		Tester for devices containing Masimo SET, 14-pin instrument connector. Tests SpO2 and PR only	No cable	Bitmos, pre-2010 RAD5/5V's